



**REPORTING NCS-PATIENT'S INTERCOM PANEL (BCP-18)**

BCP-18 is the bedside Patient's panel which enables 2-way communication between the patient and the management staff with its speaker via the microphone on the pear push switch/handset (NCB-03/05). It has 2 buttons to adjust the volume.

BCP-18 is equipped with an 8-pin military connector to ensure constant connectivity of the panel to Patient's Handset / Pear push Button despite probable jolts.

It also has a Cancel/Presence and a Nurse Call switch and provides the 12VDC to Pull Cord Switches (TB-12) and Dome Lights (OLD-19) in the room/ward through the on-board switching power supply. It can also control the patient's study light via the pear push switch.


**Technical Specifications:**

Weight	217.5 g
Dimensions	W: 227 mm H:80 mm D:40 mm
Raw Material	ABS + Polycarbonate
PowerSupply	12 V DC(via an external Power Adaptor)
Color	White
Port1	RJ45 Connectbr to Layout Ring Wiring
Port2	8-pin, Military male connector to Patient's Handset
Speaker	8 ohms, 0.5 Watt
Operating Temperature	-10 ° C ~ +50° C
Operating Humidity	30% ~ 80% RH
Mounting	Flush (Cavity Dimension: 192m x 60 mm)