

REPORTING NCS-STAFF CONSOLE (CUI-04)


CUI-04 is a desk unit with integral operating elements. It can make Direct Audio Call to each Patient Consoles (BCP-18). It can also initiate 12 emergency codes using a combination of dedicated buttons through connection to hospital Emergency /Group Summon System.

Patients can also initiate a 2-way Audio call to CUI-04 through their Service Console (BCP-18)

It has 5 indicator lamps/ Sound Alarms to notify the Nurse Call Type with a Queue management system to prioritize (Critical the highest – Emergency/Service the medium and Normal the lowest Priority) and display them in priority order on a 128 x 64 LCD display.

The history of calls and the unfinished calls can also be monitored in priority order through dedicated buttons.

Critical Calls through Critical Patient's handset/Pear Push Switch can also be set using a dedicated button on CUI-04.



Technical Specifications:

Weight	978 g		
Dimensions	W: 264 mm	H: 49mm	D: 210mm
Raw Material	Polycarbonate/ABS		
Power Supply	12 VDC		
Color	White		
Port1	DB15 Connector to Terminal Box (TDB)		
Port2	Rj45 to Hospital Emergency/Group Summon System		
Speaker	8 Ohms , 0.5 Watt		
Operating Temperature	-10°C ~ +50°C		
Operating Humidity	30% ~ 80% RH		
Mounting	Desk /Surface		