


## BASIC NCS-PATIENT PANEL (PBCN-13 & PBC-12 & PBCN-11)

PBCN-11 and PBC-12 Patient panels come with a military socket for connection to the patient pear push Nurse Call Switch (PB-03).

PBCN-11 has an additional Cancel/Presence switch whereas PBC-12 has a Nurse Call switch next to the socket.


PBCN-13 has a Nurse Call and a Cancel/Presence switch and no military socket.











## **Technical Specifications:**

| Weight                | PBCN-13: 61 g                          | PBC-12:68 g | PBCN-11:69 g |
|-----------------------|--|-------------|--------------|
| Dimensions            | W:80 mm                                | H:80 mm     | D: 20 mm     |
| Raw Material          | Polycarbonate + ABS                    |             |              |
| Color                 | White                                  |             |              |
| Connection            | 5-Pin Military Gx16 (Male)             |             |              |
| Power Supply          | 12 V                                   |             |              |
| Port                  | 5-Pin Male (Gx16-5) Military connector |             |              |
| Operating Temperature | -10°C ~ +50° C                         |             |              |
| Operating Humidity    | 30% ~ 80% RH                           |             |              |
| Mounting              | Flush/Surface                          |             |              |